

**Faculty of Biology  
and Medicine**  
Protein Analysis Facility

# Proteome quantitation by mass spectrometry

*A not-so-long, not-so-boring guide  
to data processing and interpretation*

**Unil.**

Version 5.2 – January 2026



# GOALS



## EXPLAIN

The final part of the data analysis pipeline



## DESCRIBE

The format of data



## PROVIDE

Options for downstream data processing and interpretation



## ANSWER

Frequently asked questions

*A separate guide will be available for PTM specific datasets (e.g. phosphoproteomics)*

# INDEX

4 General MS workflow

5 Things to remember

6 Different MS pipelines

7 Data analysis overview

8 Data analysis

18 Final table format

19 Unique and razor peptides

20 MaxQuant output variables

21 General facts to remember

22 Evaluation of processed output

23 Oops, it did not work


24 Publication & data submission to repositories

25 Informations resources


26 FAQ's

29 Appendix

# Proteomic workflow




# Thinks to remember *before going further*




# Different pipelines, same type of output


## Sample preparation




## MS analysis mode



## Data analysis software



## INITIAL TABLE





## FINAL TABLES & FIGURES

# Filtering to remove irrelevant data such as ...

## Reverse hits

In large datasets, a search is performed against a *decoy* database made of reverse sequences to estimate the false positive rate. Scoring parameters are then adjusted to achieve a **1% false discovery rate** (FDR, i.e. 1% decoy hits) in the final results.



## Contaminants such as keratins, trypsin, albumin

They often show up in samples either because they are part of the experimental set-up or simply common environmental contaminants.

## Proteins only identified by site

Protein groups identified **only** by a modified peptide.



## Filter on valid values


Protein groups with a high number of missing values can't be quantified with confidence – it is like trying to judge a competition when half the contestants didn't show up.



# Log2 transformation to normalize the distribution & facilitate plotting

## RAW LINEAR data

*Vast dynamic range*  
*Asymmetric distributions*



Multiplication  
Division  
0  
**FOLD change**  
 $A/B = X$

$A/B = 2$   
 $A/B = 0.5$


**Examples**

## TRANSFORMED LOG data


Sum  
Subtraction  
NaN (Not A Number)  
**log<sub>2</sub> difference**  
 $\log_2(A) - \log_2(B) = \log_2(X)$

$\log_2(A) - \log_2(B) = 1$   
 $\log_2(A) - \log_2(B) = -1$

*Manageable dynamic range*  
*More normal distributions*




# Plotting *to look at data globally*



## Medians Distributions

Detect loading  
differences

Detect measurement  
biases



## Negative & positive controls

Assess sensitivity

Biological quality control



## Correlation between replicates / Clustering

Quality control

Look for outliers,  
emerging groups


# Normalization *to compensate for global differences in sample amount*

Raw data

Assumption

Methods

Normalized data



Most of the proteome does not change (stable background)

Median subtraction is the most conservative normalization.

For IPs, data can be normalized based on intensity of a protein of interest; however, this approach may distort some datasets and lead to false positives.

Some databases, such as MaxQuant LFQ or Spectronaut Quantity values are already normalized.



# Imputation to deal with missing values

1

Quantitative proteomics data always contain **missing values**

2

Missing values could come from :

- stochasticity of data acquisition
- peptides/proteins near/below the limit of detection
- truly absent proteins

3

It is impossible to discriminate between the events in point 2.

4

## MISSING VALUES

*“0” in linear data*

*“NaN” in log data (Not A Number)*

5

Some types of data contain more missing values (LFQ, iBAQ) than others (TMT, DIA)

6

**Missing values** are a data analysis problem.

Fold changes cannot be calculated. Proteins cannot be plotted.

7

## Imputation

= substitution of missing values with values used for computation

8

Imputed values are often chosen at the low end of the distribution of the values, assuming that the reason for the missing values is very low signal intensity.

A (very) low fixed value could also be used in some cases.

9


**Imputation has to be applied with caution and on a limited scale**

Imputed data should be flagged to be easily identified (and removed) later.

# Statistical tests

*to determine whether the observed differences are statistically significant*

Example



Statistical test

Student's or Welch's tests are the most commonly used.

Statistical power

i.e. the number of replicates. This parameter is crucial, and a minimum of 3 replicates is required to ensure reliable statistical significance.

Multiple testing correction

When working with large datasets, it is important to correct p-values – for example using the Benjamini-Hochberg method – to obtain adjusted p-value and avoid falsely identifying too many results as significant.

To pass the test, the difference between groups must be larger than the variability within groups – basically, the signal has to shout louder than the noise.

# The volcano plot


a common way to represent statistical significance (y-axis) and fold-change (x-axis) in a single plot.

By using  $-\log_{10}$ , low p-values (highly significant) appear toward the top of the plot.

Significant protein groups are highlighted using a color-coded system.

- p-values < 0.05
- adj. p-values < 0.05

Arbitrary threshold for data interpretation.  
 $\log_2(\text{treatment}/\text{control}) > 1$  means treatment > 2x control in linear data



Each dot represents a protein group, and the gene names of the most relevant ones are highlighted.

# **Remove imputed values**

**The final dataset should contain only real measured values  
Imputed values are turned back into NaN (Not a Number)**

**Are the candidate proteins enriched in specific pathways, localization ... ?**

## Annotation enrichment analysis

A statistical approach is used to identify annotations that are significantly over-represented in association with a quantitative trait, typically the fold-change between 2 conditions.

### Limitations

- *Statistical bias, as different biological domains are unevenly annotated or curated*
- *Redundant or complex annotations, which can complicate interpretation*



## Functional annotation

Terms describing functional or structural features, interactions, pathways or other relevant characteristics are added to the protein groups.

Sources (not exhaustive):

<i>Gene Ontology (GO)</i>	<a href="http://geneontology.org/">http://geneontology.org/</a>
<i>KEGG</i>	<a href="https://www.genome.jp/kegg/">https://www.genome.jp/kegg/</a>
<i>Corum</i>	<a href="https://mips.helmholtz-muenchen.de/corum/">https://mips.helmholtz-muenchen.de/corum/</a>
<i>Reactome</i>	<a href="https://reactome.org/">https://reactome.org/</a>

# Table clean-up

*to get a final and comprehensive dataset*



- Columns are removed and reordered.



- You will receive a **ZIP** file containing a full reporting, plots (PNG, SVG and clickable HTML) and tab-separated tables.



- A table filtered for proteins identified by 3 or more peptides can be provided (more reliable quantitative values).

# Example of a final table, *the exact names of the columns may vary*

Quantitative values  
Could be iBAQ, LFQ,  
intensity, quantity

Includes shared peptides but assigned only to highest scoring protein group (= razor peptides).  
Peptides used for quantitation.

Difference of means of conditions.  
Equivalent to fold change, but in log<sub>2</sub> scale

Total number of peptides used for identification/quantitation.  
Includes shared peptides

Unique peptides for this protein group

Adj. p-val = p-value after multiple testing correction

Is protein group passing statistical test?  
Usually with adj p-val < 0.05

t-test p-value (uncorrected)

LOG2 scale!  
Normalized values can be negative ie lower than median

Database identifier.  
First one is the best

Use this for filtering

Total number of matched spectra.  
Can also be used for filtering/validation if column is present (SILAC or TMT)

Functional annotations  
(many terms per line)

Ctrl R1 iBAQ.log2	Ctrl R2 iBAQ.log2	Ctrl R3 iBAQ.log2	Treated R1 iBAQ.log2	Treated R2 iBAQ.log2	Treated R3 iBAQ.log2	iBAQ	Majority.protein.IDs	Protein.names	Gene.names	Fasta.headers	Peptides	Razor.unique.peptides	Unique.peptides	Ctrl R1.Razor.unique.peptides	Ctrl R2.Razor.unique.peptides	Ctrl R3.Razor.unique.peptides	Treated R1.Razor.unique.peptides	Treated R2.Razor.unique.peptides	Treated R3.Razor.unique.peptides	Mol.weight.kDa	Sequence.coverage	Score	Intensity	MS.MS.count	p.value.Treated-Control	adj. p.value.Treated-Control	log2.fold.change.Treated-Control	is.significant.Treated-Control	GOBP name	GOMF name	GOCC name	KEGG name
3.52	3.44	3.29	NaN	NaN	NaN	2.37E+08	P63328	Serine/threonine-protein p	Ppp3ca	sp P63328 PP2	22	22	14	15	14	14	0	0	0	58.6	36.3	159.18	1.34E+10	269	NaN	NaN	NaN	false	G1/; cataly nucle MAP			
5.56	5.71	5.61	-0.30	-1.12	1.66	1.66E+09	P63330	Serine/threonine-protein p	Ppp2ca	sp P63330 PP2	21	21	4	15	15	14	1	1	1	35.6	70.6	323.31	5.24E+10	580	8.3E-05	2.0E-04	-5.7	true	mitc/G-pro protei mRN			
6.22	6.54	6.03	3.39	2.17	1.50	3.22E+09	P67778;Q5SQG5	Prohibitin	Phb	sp P67778 PHI	27	27	27	14	16	16	12	7	7	29.8	73.9	323.31	1.45E+11	637	2.3E-04	4.5E-04	-3.7	true	neg; protei nucleus;n			
3.52	3.93	3.74	NaN	-1.2	NaN	3.78E+08	P67871;G3UZJ5;G3UZA	Casein kinase II subunit b	Csnk2b	sp P67871 CSK	9	9	9	8	8	8	0	1	0	24.9	59.5	153.91	1.06E+10	228	1.4E-04	3.0E-04	-4.6	true	mor chron chron Ribc			
7.53	7.80	7.75	4.5	5.41	5.23	5.53E+09	P67984	60S ribosomal protein L22	Rpl22	sp P67984 RL2	9	9	8	7	5	5	3	3	4	14.8	39.1	153.9	6.67E+10	222	2.0E-05	7.1E-05	-2.6	true	cell patter nucle Ribc			
7.46	7.62	7.49	1.96	1.35	2.18	4.53E+09	P68040	Guanine nucleotide-bindir	Gnb2l1	sp P68040 RA	26	26	26	20	18	18	6	6	7	35.1	78.5	323.31	1.50E+11	794	7.6E-08	2.8E-06	-5.8	true	regu enzym phag; Mea			
-0.28	-0.60	-1.10	NaN	NaN	NaN	2.85E+07	P68181	cAMP-dependent protein k	Prkacb	sp P68181 KA	17	4	4	2	1	1	0	0	0	40.7	31.6	12.47	1.02E+09	44	NaN	NaN	NaN	false	neu nucle nucle MAP			
7.25	7.31	7.39	2.35	2.00	2.91	3.99E+09	P68254	14-3-3 protein theta	Ywhaq	sp P68254 143	26	22	22	14	13	14	7	6	9	27.8	64.1	323.31	2.14E+11	735	3.2E-07	5.7E-06	-4.8	true	reg; bindir cyton Cell			

# Some definitions *about unique and razor peptides*

Proteins 1 and 2 have peptides in common


**Protein 1**



**Protein 2**



**DIGESTION** →




To which protein should the peptides be assigned for quantification ?


**Razor peptides**

= non-unique peptides usually assigned to the protein group with the highest number of peptides  
Protein 2 in this case

**Unique peptides**



# More definitions *MaxQuant* output variables



*Spectronaut* calculates «quantities», which are mostly similar to *MaxQuant* LFQ values

# General facts to remember



## Proteomes

**tend to change less than transcriptomes**

*The amplitudes of changes are often smaller and the changes slower.*



## Statistical power

**is crucial to find significantly different proteins**

*Power is mostly determined by the number  $N$  of biological replicates.*

*$N=3$  is the minimum for reliable (and publishable !) data.*

*$N > 3$  is recommended when quantifying mild changes, in case of high variability or for studies on individual animals.*




## Quantitation reliability

**is lower for weak signals and proteins with a low peptide count**

*These are mostly low abundance proteins.*

# Evaluation of processed output *How to define a list of candidates*



# Oupsss, it did not work out

**Sometimes, no statistically significant proteins emerge**



## Most frequent reasons:

- Small differences between conditions
- Insufficient statistical power (**N** too low)
- High variability within groups

Often more than one of these factors combined !

**It does not necessarily mean there is no difference**



It may simply indicate that the experiment was not perfect.

**How to move on ?**

- Select the best candidates (using the best fold-change & p-value, proteins identified with min. 3 peptides)
- Try to quantify in a targeted way (WB, FACS, ELISA)
- Look at annotation enrichment for global trends
- Repeat with more replicates

# Publications and data submission to repositories

Please read our guidelines on **acknowledging core facility work and authorship**.

<https://wp.unil.ch/paf/policy/>



Remember to ask for **Material & Methods** descriptions sooner rather than later.

Retrieving all relevant information about your experiment may become challenging after a few years.



Upload of raw and processed data **to repositories** is highly recommended or mandatory for most journals.

We recommend **proteomexchange.org** (part of PRIDE <https://www.ebi.ac.uk/pride/>).

We handle data preparation and upload for our users. However, we require **at least 2 weeks** prior to manuscript submission. For really complex datasets, additional time may be required. We may need to charge you for the work involved.

Submission to the data repository requires

- Raw MS data
- Project title, description
- Materials and Methods
- Protein ID and quant data

A **template for submission** to PRIDE/Proteomexchange is available on :

<https://wp.unil.ch/paf/useful-infos/>

# Ressources



LAB

**MaxQuant**

<https://maxquant.net/maxquant/>

MaxQuant and Perseus are free  
but for *Windows or Linux only!*

Tutorials can be found on YouTube

(<https://www.youtube.com/c/MaxQuantChannel/videos>)



**Perseus**

<https://maxquant.net/perseus/>



**Spectronaut**

<https://biognosys.com/software/spectronaut/>



Many **R** packages used for genomics data analysis can be applied to proteomics data after raw MS data processing (e.g. for normalization, annotation enrichment, ..).

Dedicated packages also exist.

**Why are you talking about protein groups and not proteins ?**

“Some proteins cannot be unambiguously identified by unique peptides (but have only shared peptides).

They are grouped in one protein group and quantified together, e.g. if all detected peptides of protein A also belong to protein B, A and B form one protein group.”

*Maxquant*

**Why should I do replicates ?**

The question should rather be “why did you not do replicates?”

Good science must be reproducible. Replicates provide « an estimate of between samples error». They also enhance reliability for low-abundance protein groups or subtle differences between samples.

**What is the point of doing a pilot experiment ?**

To test and optimize the different steps of both your and our experimental designs (such as proteins extraction and digestion from our side). A pilot gives an idea about feasibility and can save everyone time and money.

**Can I analyse together data from replicates that were generated and analysed at different moment ?**

Technically yes but there will be batch effects.

By far the best consistency is obtained when all samples to be compared are prepared and analysed in the same session.

**How does imputation work ?**

Please read page 12 of this manual

**Is imputation random ?**

Yes, imputation is usually done with random values taken from a low-shifted distribution, to simulate “real” data for low abundance proteins. A fixed low value can be used but this can skew the results of statistical tests, as they assume “normal” data distributions. Fixed values are used only when  $N=1$ .

**Is imputation reproducible?**

Not 100%, since the values used are random-generated. If the imputation is repeated, the imputed values will be slightly different.

**Imputation is random and not reproducible :  
how is it reliable ?**

Since the values are all low and in a small range, the impact of the randomness of the imputed data on the final results should be limited. We also do not perform calculations on totally imputed data, i.e. when comparing two groups, at least in one group there should be real values.

**How do I identify outlier samples ?**

The PCA (principal component analysis) plot is the best tool to spot outliers. It also enables us to evaluate sample homogeneity within groups and differences between conditions.

**Why are data filtered using peptides number ?**

We assume that protein groups identified with few peptides may represent low-abundance proteins / may be less reliable.

We typically prioritize protein groups identified with minimum 2 Razor.unique peptides (i.e.  $>1$  razor.unique peptides).

**Why data are filtered using a certain number of valid values into a group (when replicates are present)?**

We assume that protein groups detected in only some replicates of a group may represent low-abundance proteins or may be less reliable. We typically prioritize protein groups consistently identified across all replicates within at least one sample group.  
*For instance, when working with triplicates, we retain protein groups that have at least 3 valid intensity values within a given group.*

**What is the minimum (fold) change that I can consider as "real" ?**

As a rule of thumb, fold changes smaller than  $\pm 0.5$  in  $\log_2$  should not be considered. A safe minimal threshold is  $\pm 1 \log_2$ , a two-fold change in linear data. Annotation enrichment analysis is a way to highlight global changes that are of low magnitude.

# **APPENDIX**

***Additional information about annotation enrichment***

# Annotation enrichment *Basic concepts*

## Annotation

Terms are added to each protein in the list of interest.

*ex: GO terms*

## Sorting

The list is sorted based on a relevant experimental parameter.

*ex: fold change between the conditions*


## Statistics

An algorithm is applied to assess whether each annotation term is significantly enriched at either the top or bottom of the list, indicating a tendency to be associated with proteins that either increase or decrease between conditions.

## Correction

A p-value is calculated and should be corrected for multiple testing in case of large datasets with many annotation terms.

Enrichment plot: GO:0006457



The algorithm used in our analyses is described in Cox & Mann (2012) *BMC Bioinformatics*, 13 (Suppl 16), S12. <https://doi.org/10.1186/1471-2105-13-S16-S12>. It is conceptually similar to the widely used GSEA algorithm.

# Annotation enrichment *Result table*

Parameter used for sorting the table.  
Usually  $\log_2$  (fold change)

P-value for enrichment of this term

Number of protein groups in the list with this annotation type

Corrected p-value (Benjamini-Hochberg method).  
A cut-off of 0.02 is usually applied to generate this list.

Column	Type	Name	Size	Score	P value	Benj. Hoch. FDR	Mean	Median
log2.fold.change.Treated_Control	GOBP name	drug transmembrane tr	6	0.92	1.01E-04	0.005	3.315	3.83
log2.fold.change.Treated_Control	GOMF name	polyamine transmemb	8	0.74	2.71E-04	0.012	3.284	3.66
log2.fold.change.Treated_Control	GOBP name	drug transport	7	0.89	4.68E-05	0.003	2.997	3.60
log2.fold.change.Treated_Control	GOMF name	drug transmembrane tr	7	0.89	4.68E-05	0.002	2.997	3.60
log2.fold.change.Treated_Control	GOBP name	mannose transport	7	0.96	1.02E-05	0.001	3.781	3.50
log2.fold.change.Treated_Control	GOMF name	fructose transmembran	7	0.96	1.02E-05	0.001	3.781	3.50
log2.fold.change.Treated_Control	GOMF name	iron ion transmembran	6	0.84	3.51E-04	0.015	2.313	2.07
log2.fold.change.Treated_Control	GOCC name	intrinsic to plasma me	24	0.56	2.68E-06	0.000	1.971	1.97
log2.fold.change.Treated_Control	GOCC name	anchored to membrane	10	0.77	2.86E-05	0.001	1.824	1.97
log2.fold.change.Treated_Control	GOBP name	copper ion transport	10	0.74	5.34E-05	0.003	2.085	1.95
log2.fold.change.Treated_Control	GOCC name	integral to plasma men	21	0.58	4.43E-06	0.000	2.071	1.91
log2.fold.change.Treated_Control	GOCC name	vacuolar lumen	6	0.84	3.97E-04	0.007	2.646	1.70
log2.fold.change.Treated_Control	GOMF name	cofactor transporter ac	8	0.77	1.64E-04	0.008	1.926	1.65
log2.fold.change.Treated_Control	GOMF name	amine transmembrane	32	0.49	2.23E-06	0.000	1.522	1.56
log2.fold.change.Treated_Control	GOBP name	cofactor transport	10	0.67	2.61E-04	0.010	1.542	1.52

GOBP : Biological process  
GOCC: Cellular Compartment  
GOMF: Molecular fonction

Terms are often redundant, due to hierarchical organization of GO terms

Mean resp. Median of «T-test differences» for each GO term, calculated on the protein groups considered (Size).  
In  $\log_2$ ,  $>0$  means increasing in Treated and *vice versa*.

Can be seen as normalized version of mean/medians.  
Values are always between -1 and +1, centered at 0.

# Annotation enrichment

*How to retrieve proteins associated with a specific GO annotation ?*


Column	Type	Name	Size	Score	P value	Benj. Hoch. FDR	Mean	Median
log2.fold.change.Treated_Control	GOBP name	drug transmembrane tr	6	0.92	1.01E-04	0.005	3.315	3.83
log2.fold.change.Treated_Control	GOMF name	polyamine transmemb	8	0.74	2.71E-04	0.012	3.284	3.66
log2.fold.change.Treated_Control	GOBP name	drug transport	7	0.89	4.68E-05	0.003	2.997	3.60
log2.fold.change.Treated_Control	GOMF name	drug transmembrane tr	7	0.89	4.68E-05	0.002	2.997	3.60
log2.fold.change.Treated_Control	GOBP name	mannose transport	7	0.96	1.02E-05	0.001	3.781	3.50
log2.fold.change.Treated_Control	GOMF name	fructose transmembran	7	0.96	1.02E-05	0.001	3.781	3.50
log2.fold.change.Treated_Control	GOMF name	iron ion transmembran	6	0.84	3.51E-04	0.015	2.313	2.07
log2.fold.change.Treated_Control	GOCC name	intrinsic to plasma me	24	0.56	2.68E-06	0.000	1.971	1.97
log2.fold.change.Treated_Control	GOCC name	anchored to membrane	10	0.77	2.86E-05	0.001	1.824	1.97
log2.fold.change.Treated_Control	GOBP name	copper ion transport	10	0.74	5.34E-05	0.003	2.085	1.95
log2.fold.change.Treated_Control	GOCC name	integral to plasma men	24	0.58	4.43E-06	0.000	2.071	1.91
log2.fold.change.Treated_Control	GOCC name	vacuolar lumen	6	0.84	3.97E-04	0.007	2.646	1.70
log2.fold.change.Treated_Control	GOMF name	cofactor transporter ac	8	0.77	1.64E-04	0.008	1.926	1.65
log2.fold.change.Treated_Control	GOMF name	amine transmembrane	32	0.49	2.23E-06	0.000	1.522	1.56
log2.fold.change.Treated_Control	GOBP name	cofactor transport	10	0.67	2.61E-04	0.010	1.542	1.52
log2.fold.change.Treated_Control	GOCC name	fungus-type vacuole lur	5	0.83	0.00133107	0.020	2.784	1.45
log2.fold.change.Treated_Control	GOBP name	transition metal ion tra	20	0.51	1.85E-06	0.000	1.359	1.40
log2.fold.change.Treated_Control	GOBP name	iron ion transport	11	0.62	4.15E-04	0.014	1.577	1.40
log2.fold.change.Treated_Control	GOCC name	organellar large riboso	41	0.67	1.16E-13	0.000	1.168	1.33
log2.fold.change.Treated_Control	GOCC name	mitochondrial large rib	41	0.67	1.16E-13	0.000	1.168	1.33
log2.fold.change.Treated_Control	GOMF name	amino acid transmemt	27	0.40	3.62E-04	0.015	1.107	1.29

How to get these 24 proteins ?

# Annotation enrichment

*How to retrieve proteins associated with a specific GO annotation ?*

Open the table in Excel



Select the column that contains the annotation of interest  
Ex: GOBP



Go to «data» and click «Filter»


The screenshot shows the Microsoft Excel interface. The 'DATA' tab is active, and the 'Filter' button in the ribbon is highlighted with a red box. Below the ribbon, the formula bar shows 'GOBP name'. The data table below has the following columns: A (Ratio H/L normalized), B (Ratio H/L normalized), C (Amino acid), D (Charge), E (T-test Significant t p005), F (GOBP name), G (GOMF name), H (GOCC name), I (id), J (Localizati on prob), K (PEP), L (Score), and M (Delta score). The 'GOBP name' column is highlighted in green, and a dropdown arrow is visible in cell F1. The table contains 15 rows of data, with rows 56-70 visible.

	A	B	C	D	E	F	G	H	I	J	K	L	M
	Ratio H/L normalized	Ratio H/L normalized	Amino acid	Charge	T-test Significant t p005	GOBP name	GOMF name	GOCC name	id	Localizati on prob	PEP	Score	Delta score
1	9889	9890											
56	0.498	0.693	Y	2		biological	regulation cell part;c		45	1	2.88E-08	121.83	109.25
57	0.401	0.750	Y	2		apoptosis	enzyme ac cell juncti		317	1	0.02792	43.083	15.558
58	0.190	0.954	Y	3		biosynthe	binding;ca cell part;c		424	0.99989	3.97E-07	71.692	67.747
59	0.225	0.915	Y	3		assembly	binding;rik Cajal body		383	0.99518	0.00031	72.652	72.652
60	0.313	0.824	Y	2		alcohol m	catalytic a brush bord		354	1	0.013	42.001	42.001
61	-0.683	1.809	Y	4		axon guid	binding;ca adherens		474	1	2.73E-15	87.659	87.58
62	-0.013	1.113	Y	2		biological	binding;ca adherens		94	0.99971	0.00587	61.999	18.753
63	0.329	0.752	Y	2		apoptosis	enzyme ac cell juncti		319	1	5.22E-09	92.063	85.024
64	-0.121	1.191	Y	3		cell cycle	cytokinesis cell divisio		332	1	4.97E-08	87.18	87.18
65	-0.210	1.213	Y	3					338	0.99245	7.02E-06	49.089	49.089
66	-1.370	2.332	Y	2		antigen p	binding;ca cell part;c		305	1	0.00111	86.794	72.637
67	-0.079	1.029	Y	3		biosynthe	binding;ca cell part;c		423	1	7.16E-05	62.924	56.761
68	-0.143	1.084	Y	3		actin cyto	adenyl nu cell juncti		109	1	2.39E-20	94.384	90.986
69	-0.131	1.043	Y	3			receptor a cell part;c		397	1	2.53E-21	102.6	102.6
70	-0.043	0.947	Y	2		alcohol bi	binding;ca cell part;c		112	1	3.17E-05	96.331	91.88


# Annotation enrichment

How to retrieve proteins associated with a specific GO annotation ?

Go to the pull-down menu



Select «text filters» and «contains»



Enter the name of the annotation you are searching for in the «contains» field

The screenshot shows an Excel spreadsheet with a data table. The columns are labeled as follows: A: Ratio H/L normalized 9889; B: Ratio H/L normalized 9890; C: Amino acid; D: Charge; E: T-test Significance p005; F: GOBP name; G: GOMF name; H: GOCC name; I: id; J: Localization prob; K: PEP; L: Score; M: Delta score; N: Score for localization; O: Mass error [ppm]; P: Intensity L; Q: Intensity H; R: Intensity H; S: Localization on prob 9889; T: Localization on prob 9890; U: Ratio H/L count 9889; V: Ratio H/L count 9890; W: Position. A 'Text Filters' menu is open over the 'GOBP name' column, showing options like 'Sort A to Z', 'Sort Z to A', 'Sort by Color', 'Clear Filter From "GOBP name"', 'Filter by Color', and 'Text Filters'. The 'Text Filters' menu is further expanded to show 'Contains...', 'Does Not Contain...', and 'Custom Filter...'. A 'Custom AutoFilter' dialog box is open in the foreground, showing the 'Show rows where:' field with 'GOBP name' selected, and the 'contains' operator chosen. The 'And' radio button is selected. The 'OK' and 'Cancel' buttons are visible at the bottom of the dialog box.



**Faculty of Biology  
and Medicine**  
Protein Analysis Facility

**That's it !**

*Hope it was useful !*

**Unil.**